sōna

AGM LFA

Aspergillus Galactomannan Lateral Flow Assay

Aspergillus Galactomannan

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Asp LFA Publication Findings

The clinical performance of the IMMY sona Aspergillus Lateral Flow Assay in serum samples from cases classified as proven/ probable/ chronic IA/IFD (n=32) using EORTC/MSG criteria and control patients with no evidence of IFD (n=100).

Sensitivity	Specificity	PPV	NPV
97%	98%	94%	99%

Lewis W., et al. An evaluation of the performance of the IMMY sōna Aspergillus GM LFA when testing serum to aid in the diagnosis of IA. J Clin Micro. 2020; 00053-20. In serum samples, the AGM LFA will correctly generate a positive result for 97% of patients with IA, and a negative result for 98% of patients without IA. The PPV shows that among patients with a positive AGM LFA result, the probability that they truly have IA is 94%. NPV shows that among patients with a negative AGM LFA result, the probability that they do not have IA is 99%.

Asp LFA Publication Findings

The diagnostic performance of the IMMY sona Aspergillus Lateral Flow Assay in BAL samples from hematology patients classified as proven IPA (n=11) using EORTC/MSG criteria and control patients with no evidence of IPA (n=117).

Sensitivity	Specificity	NPV
91%	92%	99%

Mercier T, et al. Lateral flow assays for diagnosing invasive pulmonary aspergillosis in adult hematology patients: A comparative multicenter study. Med Mycol. 2019. In BAL samples, the AGM LFA will correctly generate a positive result for 91% of patients with IA, and a negative result for 92% of patients without IA. The NPV shows that among patients with a negative AGM LFA result, the probability that they do not have IA is 99%.

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simple Procedure

SPECIMEN PREPARATION

- 1 Obtain 2 test tubes for each specimen.
 - 1 screw cap, heat resistant centrifuge tube for the dilution
 - 1 flat-bottom tube for running the test
- 2
 - Transfer 300 µl specimen in the screw cap, heat resistant centrifuge tube 1
 - Add 100 µl Sample Pretreatment Buffer to tube 1 (*(vortex as needed)*)

- 3
 - Place tube 1 on heat block for 6-8 minutes at 120 C
- 4
 - Centrifuge tube 1 at 10,000 - 14,000 x g for 5 minutes

- 5
 - Transfer 80 µL from tube 1 to tube 2
- 6
 - Add 40 µL of Aspergillus GM LFA Running Buffer to tube 2

- 7
 - Insert strip (down) Incubate for 30 min.
- 8
 - Read Test
 - 1 line = negative
 - 2 lines = positive

RUN TEST

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Check Diagnostics GmbH

Bergkoppel 11 23843 Bad Oldesloe

Tel.: 04531 42835 00

Fax.: 04531 42835 10

check@check-dx.de

www.check-dx.de