

Enviroscreen

PreMoist Swab Kits



- > Neutralise disinfectant residues
- > Integral solution for HACCP
- > Entire sample can be analysed
- > Effective hygiene monitoring solution
- > Ideal for validation of TACCP
- > Sample integrity & security
- > No loose packaging in critical areas
- > Choice of swab materials
- > Strengthening brand protection

No peel pouch & ready to use for convenience
Simple & effective sampling kit
Ideal for environmental and hygiene monitoring of
equipment including critical surfaces

ISO 9001:2008
ISO 18593 compliant

12 months shelf life
Sterilised by
Gamma Irradiation



Technical Service
Consultants Ltd

Technical Service Consultants Ltd

Microbiology House, Fir Street, Heywood, Lancashire OL10 1NW, UK
T : +44 (0)1706 620600 F : +44 (0)1706 620445
E : sales@tscswabs.co.uk W : www.tscswabs.co.uk

PreMoist Swab Sampling Kits

PreMoist sampling swabs offer end-users a convenient and quick sampling method to confirm the presence of bacteria with no additional preparation required. PreMoist Swabs can be used for environmental and hygiene monitoring, food security and cleaning validation due to the inclusion of neutralisers.

In many cases a surface will contain traces of disinfectants and/or sanitisers. Left alone and over time, these cleaning fluid traces will not kill the microorganisms, but will just prevent any growth (bacteriostatic effect). TSC Neutralising Buffer will neutralise most common disinfectants including peroxides and derivatives, phenols, quaternary ammonium compounds (QAC's), amphoteric, chlorines and associated halogens, providing end-users with total testing assurance.

Surface Sampling Swabs

- > Easy direct sampling
- > Choice of viscose and foam tips on white or blue polystyrene shafts
- > Foam Tip Swab increases sample collection & release
- > Define sampling areas by using a sterile sampling template (TS/15-T)
- > Labelled tubes for easier sample processing and identification
- > Tamper evident seals & shatterproof tubes offer total product confidence









Download Certificates of Conformity from
www.tscswabs.co.uk

How to use guides available:
www.tscswabs.co.uk/sampling

Sampling Areas

- > Entire sample can be analysed
- > Ideal for sampling a wide range of environments and critical surfaces > Exposes potential microorganism contamination
- > Useful in hygiene surveillance (HACCP) and environmental monitoring > Food security / TACCP validation
- > Protecting your brand and preventing costly recalls



Code	Description	Qty	Recommended Use		
			Environmental	Food & Beverage	Clinical & Veterinary Infection Control
TS/6-61	 <ul style="list-style-type: none"> > PreMoist viscose tip swab in TSC Neutralising Buffer > Will neutralise sanitisers and disinfectants > White polystyrene shaft > Suitable for most surfaces 	160	✓	✓	✓
TS/6-61BP	 <ul style="list-style-type: none"> > PreMoist viscose tip swab in TSC Neutralising Buffer > Blue polystyrene breakpoint shaft for high visibility > Sample easily transferred to smaller vials using breakpoint > Will neutralise sanitisers and disinfectants > Suitable for critical surfaces in a production area 	160	✓	✓	
TS/6-62	 <ul style="list-style-type: none"> > PreMoist viscose tip swab in deionised water > White polystyrene shaft > Suitable for when Neutralising Buffers are not appropriate > Ideal for sampling sensitive equipment 	160	✓		✓
TS/6-63BP	 <ul style="list-style-type: none"> > PreMoist foam tip swab in TSC Neutralising Buffer > Foam tip swab offers enhanced collection and release of sample > Will neutralise sanitisers and disinfectants > White polystyrene breakpoint shaft > Sample easily transferred to smaller vials > Suitable for most production surfaces 	160	✓	✓	✓

Im Vertrieb von:

Check Diagnostics GmbH

Bergkoppel 11 23843 Bad Oldesloe

Tel.: +49 4531 42835 00 Fax.: +49 4531 42835 10

check@check-dx.de

www.check-dx.de

